

GradyVet: DUS - 4000

The Future of Veterinary Ultrasound



- Ultra-Compact
- Full digital beam-former
- Harmonic imaging (THI)
- Built-in lithium battery (re-chargeable by car)
- Water-proof keyboard cover (option)
- i-Image
- SRA
- Chroma, Trapezoid imaging, Compound image



Canine Liver



Canine Intestine



Bovine Gender



Bovine Follicle



FEATURES



www.GradyMedical.com
Phone: 1.800.800.2585
Email: sales@gradymedical.com


DUS-4000 Probes



D7L40L - 7.5MHz Linear
Applications: Small animal ABD, OB, GYN, Cardiac




D3C60L - 3.5MHz Convex
Applications: Small animal ABD, OB, GYN, Cardiac




D6C15L - 6.0MHz Micro-Convex
Applications: Small animal ABD, OB, GYN, Cardiac



D5C20L - 5.0MHz Micro-Convex
Applications: Small animal ABD, OB, GYN, Cardiac



Compound Imaging



Before Compound imaging



After Compound imaging

The traditional Compound imaging technologies improve the image quality, but lose the frame rate, which affect diagnosis convenience and efficiency. Our Compound imaging has integrated multiple key factors like frequency, focus and time, but will not sacrifice the frame rate.



Various measurements



Detailed reports



EasyView image archive system

WORKFLOW

The multi-purpose user presets, comprehensive measurement & report system, built-in EasyView image archive system, quick image storage/retrieve/transfer, one-button direct print, makes the complete workflow better than what you can dream of.